

ZXQ-1 Automatic Mounting Press Metallographic



PRODUCT DESCRIPTION

ZXQ-1 automatic metallurgical sample mounting press adopts automatic mechanical pressure. The mounting press is controlled by the program, which realizes the one-key operation of heating , loading, holding, cooling and unloading process. The sample mounting is automatically completed with sound and light alarm function, without the need of the operator on duty.



PRODUCT PARAMETER

Model	ZXQ-1
Method	Automatic
Mould Diameter	Φ25mm or Φ30mm or Φ40mm

Cavity Height	115mm, effective 60mm
Cooling Method	External water source cooling
Heating Power	Φ25/Φ30mm 1.2Kw Φ40mm 1.6Kw
Heating Temperature	50-180℃ (The default is 120℃)
Heating Time	1-7min(The default is 4min 30s)
Cooling Time	1-30min(The default is 2min 30s)
Power Supply	220V,50Hz, 10A

COMPANY INFORMATION



Chongqing Scope Instrument Co., Ltd. is a Chinese leading supplier who focus on laboratory instrument, and our team has 20 years rich experience in the field of laboratory instrument. We have perfect relationship with most of leading manufacturers in the field of laboratory instrument in China, and supply perfect service and different kinds of quality laboratory instrument from different leading manufacturers to our customers. We are now serving for more than 50 customers from Europe, Asia, America, Africa and Australia, and supply high quality laboratory instrument and perfect service to them. Also, we established good relationship with more than 150 customers in worldwide market. Our laboratory instrument are used widely in the field of education, medicine, petroleum, geological mineral and chemical industry and electronic industry.



